



# STRETCHING YOUR CLINICAL OPTIONS WITH "INFINITI"

OUR INFINITI LINE OF PRODUCTS PROVIDES UP TO 10 CHANNELS OF WIRELESS BIOFEEDBACK






## BIOFEEDBACK CAN MAKE A WORLD OF DIFFERENCE FOR REHAB

Biofeedback keeps patients motivated using the latest rehab techniques from graphical and audio feedback, to automated training and review protocols.

### ONE CLICK START UP PROVIDES:

- Surface Electromyography Biofeedback
- Range of Motion Training
- Respiration Training
- Heart Rate Variability Training
- Skin Conductance and Peripheral Temperature Biofeedback



## COMPREHENSIVE VIEW OF PATIENT REHABILITATION

Choose from a wide variety of assessment protocols with specialized reports:

- Dynamic and Static SEMG Assessment
- Dynamic and Static Range of Motion Assessment (with optional SEMG)
- Muscle Testing (with optional SEMG)
- Algometry
- Video recording and review to help patients follow protocols
- Automated protocols guide treatments

## TELE-INFINITI CF BLUETOOTH CARD

Roam freely and send signals up to 300 feet (100 meters)




## MYOscan-Z SENSOR - ASSURES GOOD CONTACT

A unique EMG sensor with built-in impedance checking for high signal quality

## INTERFACE TO ISOKINETIC MACHINES

Monitor EMG and other signals from Biodex and other isokinetic machines



**Thought Technology Ltd.**

<http://www.thoughttechnology.com>

Tel: (800) 361-3651 or (514) 489-8251

**MOTIVATION - COMPLIANCE - SUCCESS!™**

MAR911-00