

EVERYTHING YOU NEED TO DO TWO OR FOUR CHANNEL Z-SCORE BIOFEEDBACK:

- ProComp2, ProComp5, ProComp Infiniti or FlexComp Infiniti with BioGraph Infiniti 6.0
- Z Score 6 (SA7962)
- Z Score add-on CD with Dr. Thatcher's normative database (SA7908)
- 2 or 4 EEG-Z sensors (T9305Z)
- A 2-channel (T8760) or 4-channel (T8761) connectivity kit

The 4 Channel Connectivity Kit:

- Our most flexible TT-EEG kit. Used in conjunction with four TT-EEG Monopolar/Bipolar Electrode kits
- Connect 1, 2, 3, or 4 monopolar or bipolar channels, with the same reference for all channels.
- Choose between linked ear or single ear reference, and between ear clip or cup electrode active ground.

Please note that the 4 Channel Connectivity Kit can be used in conjunction with up to 4 EEG-Z sensors and 4 monopolar/bipolar kits (T8750). Each EEG-Z sensor includes one T8750 kit.

BioGraph® Infiniti Platform

- **Stability.** BioGraph® Infiniti is the most reliable and accurate software in the biofeedback market. Never lose data again, with our exclusive data lock technology and pause on disconnect features.
- **Flexibility.** A simple right click allows you to change settings on the fly - such as bandwidths, sound feedback and animations. Observe changes in other physiology, using any of our other sensors by customizing your screens with our FREE Developer Tools.



Thought Technology Ltd.
<http://www.thoughttechnology.com>
Tel: (800) 361-3651 or (514) 489-8251
Fax: (514) 489-8255




Strict adherence to the
IEEE Recommendations for
Neurofeedback Systems

Join us on...



Z-Score⁶ SUITE



Enhance Your Z-Score Biofeedback

SIMPLE Z-SCORE ASSESSMENT & TRAINING

EVERYTHING YOU NEED FOR 2 AND 4-CHANNEL Z-SCORE BIOFEEDBACK

THOUGHT TECHNOLOGY BRINGS YOU THE MOST EASY TO USE, POWERFUL, AND FLEXIBLE Z-SCORE BIOFEEDBACK SYSTEM ON THE MARKET.

- 1. NO MORE GUESSING.** Train using simple metrics that report on multiple z-score values simultaneously. Choose between 4 preset indices, modify them easily or make your own.
- 2. SAVE TIME, AVOID THE HEADACHE.** Streamlined performance and flexibility for all your 2 and 4-channel Z-Score Biofeedback activities.
- 3. KEEP EVERYTHING IN ONE PLACE.** Store and Load your assessment screens, training screens, statistics, and review screens with our exclusive QuickStart one-button session control.



Who Knew Z-Score Biofeedback Could Be This Much Fun:

- Bright, colorful screens.
- Single and dual monitor displays, DVD, animations & sound feedback.
- Powerful real-time trending of means, both in-session and post-session.
- Seamless changes of graphical settings from built-in mean right on the screen.
- Reviewing and reporting on screens that include efficient graphical representations of session trends.

Why Z-Score 6?

NEW! Automatic Real-Time Artifact Rejection

- Calibrate real time artifact rejection for each individual client.
- Store client-specific data so it only has to be done once, or recalibrate every time.
- Control artifact source by choosing any combination of electrodes, and easily turn the feature on or off.

Artifact settings apply to all EEG data simultaneously, either feedback (in real time) or statistics (offline review).

NEW! Layout and design

Access all individual z-scores in the most user-friendly interface on the market.

NEW! All-Z-score screen instrument with real time colour change.

NEW! Helpful sample sessions included to demonstrate the feedback and artifact rejection functions.

