

YOU ASKED
WE DELIVERED



BioGraph[®]
INFINITI
VERSION 6.1

NEW FEATURES AT A GLANCE

- Redesigned DVD Instrument -
- HTML and PDF Report Formats -
- Automatic Session Pause -
- Mute Function on Individual Instruments -
- Many Overall Improvements -

BIOGRAPH INFINITI

REDESIGNED DVD INSTRUMENT

The DVD feedback instrument was rebuilt from scratch, in order to make it completely independent of third-party audio video encoding tools (codecs) and improve its reliability and stability. Two new feedback actions were also added for "Brightness" and "Saturation".

HTML AND PDF REPORT FORMATS

Previous versions of BioGraph required Microsoft Word to generate session and trend reports. You can now generate reports in HTML – viewed in your default web browser - or in Portable Document Format (PDF) – viewed with your default PDF reader.

AUTOMATIC SESSION PAUSE

A new option to automatically pause the session after a set period of time was added to increase the flexibility of open display sessions. With this global setting, you can easily define to run trial-based sessions with breaks between pre-set periods of training.

MUTE FUNCTION ON INDIVIDUAL INSTRUMENTS

Now, feedback sounds on individual instruments can be turned ON or OFF with a simple right click selection. Now, you can rapidly and easily personalise the audio feedback to your client's needs without interrupting the training.



OVERALL IMPROVEMENTS

New main screen menu items provide one – click access to online training videos and software updates. This enhances our user's out- of the box experience.

WHAT ARE YOU WAITING FOR?

We delivered our most powerful and user friendly BioGraph yet!
What are you waiting for to upgrade?



www.thoughttechnology.com
TEL: (800) 361-3651 | FAX: (514)489-8255



MAR1095-00