



# Z-Score Biofeedback Made Easy

Thought Technology's products for z-score assessment and training combine to create the most easy to use, powerful, and flexible Z-Score Biofeedback system on the market.

Try it today, and find out why BioGraph® Infiniti, now with Z-Score Biofeedback Technology, is the most widely used biofeedback software in the world.

# Simple Z-Score Assessment & Training

*Everything you need for 2 and 4-channel Z-Score Biofeedback*

- **No more guessing.** Train using simple indices that report on multiple z-score values simultaneously. Use any of 4 preset indices, modify them easily or make your own.
- **Measure standard bands of EEG activity at the same time.**



## Quick Start Icon


Try the new **one-button session loading** that includes everything you need:

- Assessment screens, training screens, statistics, and review screens.
  - Access all this directly via a single desktop icon!
- 
- Use **bright, colorful screens** that include popular BioGraph<sup>®</sup> Infiniti features such as single and dual monitor displays, DVD, animations & sound feedback, and powerful real-time trending of means in-session and post-session.



## Easy Access to Settings

Enjoy easy access to simple Graphical User Interfaces. Access graphical settings via right click menus directly on the biofeedback screen.



## Report Screens

Review and report using screens that include efficient graphical representation of session trends.



## BioGraph® Infiniti Platform

- **Stability.** BioGraph® Infiniti is the most reliable and accurate software in the biofeedback market. Never lose data again, with our exclusive crash-proof technology and pause-on disconnect features.
- **Flexibility.** A simple right click allows you to change settings on the fly - such as bandwidths, sound feedback and animations. Observe changes in other physiology, using any of our other sensors by customizing your screens with our FREE Developer Tools.

# Electrode kits tailored for Z-Score Biofeedback ... and ready for all your EEG applications

Thought Technology's new line of EEG electrodes delivers high quality performance, ease-of-use and maximum flexibility. Use efficient 2 and 4-channel kits designed for **Z-Score Biofeedback**.



- All electrodes are available in **gold** or **sintered silver-silver/chloride\***. Several colors available!
- Ready-to-use kits have been designed with your applications in mind and allow for multiple connection configurations
- **Individual cables** are also available to create any configuration you like
- All cables feature **standard DIN connectors** for compatibility with Infiniti and **other EEG systems**
- Kits include new **reference ear clip**.
  - \* Comfortable, sturdy easy to apply and quick to clean. Includes standard DIN connector for use with any EEG equipment.
  - \* Also available separately as a useful add-on to third party **electrode caps!**

## Two and four-channel connectivity kits



- Apply fewer electrodes and record more accurately
- Use single or **linked ears reference**, common ground and multiple configurations.
- All with the same cable set.
- Ideal for **Z-Score Biofeedback!**

\* Required for all DC recordings to minimize electrode offset and drift.



**Thought Technology Ltd.**

Tel: 1-800-361-3651 • 514-489-8251 Fax: 514-489-8255  
<http://www.thoughttechnology.com>

MAR937-00