



Thought Technology Ltd.

2180 Belgrave Avenue, Montreal Quebec, H4A 2L8 Canada
 Tel: 800-361-3651 (514) 489-8251 Fax: (514) 489-8255
 Email: mail@thoughttechnology.com
www.thoughttechnology.com

BioGraph Infiniti Computer Requirements

Recommended

- IBM PC compatible, AMD Athlon XP 3000 or higher, Pentium P4 CPU speed 3 GHz or higher or equivalent mobile Laptop CPU.
- Desktop or Laptop with two monitor capability
- Windows 2000/XP Professional or Home edition.
- 50 - 60 gigabytes hard disk space for video recording and processing.
(The software needs 2.5 gigabytes to install and run of available hard drive space)
- Memory, 512 MB of RAM or more
- CD ROM or DVD drive
- SVGA graphic card (1024 x 768) or higher resolution adapter & monitor
- 32 bit compatible sound card & speakers
- 1 to 4 USB ports, depending on the desired number of Pro/Flex Infiniti encoders
- 1 to 4 serial com ports, depending on the desired number of ProComp+ encoders.
- Mouse or compatible pointing device
- MS Word 97 or higher (for printing purposes)
- Compact Flash Reader (For use with compact flash card only)
- Webcam 30 frames per second (for video purposes only)

NOTE: When using certain more complex screens, such as Coherence, you must adhere to the Recommended Computer Requirements.

CardioPro 2.0 Computer Requirements

Recommended	Minimum
<ul style="list-style-type: none"> • Computer laptop or desktop • Windows '98, 2000, XP • Processor, Pentium 4 2.4 GHz • Memory, 512 MB of RAM • Hard Drive, 40 GB • CD ROM • 15" SVGA monitor, Video card with 32 or 64 MB • Sound Blaster AWE 32 sound card & speakers • 1 free serial port, COM1 or COM2 • Mouse or compatible pointing device 	<ul style="list-style-type: none"> • Computer laptop or desktop • Windows '98, 2000, XP • Processor, Pentium II 350MHz • Memory, 128 MB of RAM • Hard drive, 5 GB • CD ROM • 15" SVGA monitor, video card with 4 MB • MB • Sound Blaster compatible sound card & speakers • 1 free serial port, COM1 or COM2 • Mouse or compatible pointing device



Thought Technology Ltd.

2180 Belgrave Avenue, Montreal Quebec, H4A 2L8 Canada
Tel: 800-361-3651 (514) 489-8251 Fax: (514) 489-8255
Email: mail@thoughttechnology.com
www.thoughttechnology.com

Continance Software Computer Requirements (for MyoTrac3 and MyoTrac3G devices)

Recommended	Minimum
<ul style="list-style-type: none">• Computer Desktop• Windows'98,2000• Processor, Pentium 4 2.4 GHz• Memory, 512 MB of RAM• Hard Drive, 40GB available disk space• Floppy Drive• 17" SVGA monitor,• Video card with 32 or 64 MB• Sound Blaster AWE 16bit sound card & speakers• 1 free serial port, COM1 or COM2• Mouse or compatible pointing device	<ul style="list-style-type: none">• Computer laptop or desktop• Windows 95, 98,2000• Processor, P2 350• Memory, 128 MB of RAM• Hard drive, 5MB available disk space• Floppy Drive• 15" SVGA monitor,• Video card with 1 MB• Sound Blaster compatible sound card & speakers• 1 free serial port, COM1 or COM2• Mouse or compatible pointing device