



FlexComp Infiniti



ProComp Infiniti



ProComp5



ProComp2



	FlexComp Infiniti	ProComp Infiniti	ProComp5	ProComp2
Product Number	T7555M	T7500M	T7525M	T7400M
Channels	10	8	5	2
Compatible Sensors	EEG, SCP/EP, EMG, EKG, BVP(Blood Volume Pulse), Respiration, Temperature, Skin Conductance, Inclinometer, Algometer/Force, Bend, Foot Switch, Hand Push Button, Voltage Isolator, Force Adapter, Goniometer/ Torsiometer Adapter, HEG(Hemoencephalography)	The same as FlexComp Infiniti	The same as FlexComp Infiniti	The same as FlexComp Infiniti
Sensor type	External sensors	External sensors	External sensors	1 Internal EEG on Channel A External sensors on B, C & D
Impedance check (with EEG-Z sensor)	All channels	All channels	All channels	Channel B only
Sampling Rate	2048 samples/sec on all 10 channels	2048 samples/sec on Channel A and B, 256 samples/sec on Channel C, D, E, F, G, H	2048 samples/sec on Channel A and B, 256 samples/sec on Channel C, D, E	256 samples/sec on Channel A and B, 32 samples/sec on Channel C and D
Connection to PC	Fiber Optic to USB through TT-USB interface unit Or Compact Flash Module Bluetooth Tele-Infinity	The same as FlexComp Infiniti	The same as FlexComp Infiniti	Fiber Optic to USB (TT-USB) or COM port (Pro-SB) or Tele-Infinity
Software and Suites	BioGraph Infiniti, Developer Tools, EEG Suite, Physiology Suite, and DynaMap Suite	BioGraph Infiniti, Developer Tools, either EEG Suite or Physiology Suite,	The same as ProComp Infiniti	The same as ProComp Infiniti
Weight	200g (7oz)	200g (7oz)	200g (7oz)	40 g (1.4oz)
Size	5.1" x 3.7" x 1.5" (130 mm x 95 mm x 37 mm)	The same as FlexComp Infiniti	The same as FlexComp Infiniti	2.2" x 2.8" x 0.75" (55 mm x 71 mm x 19 mm)
Self-calibration	Yes	Yes	Yes	No
Compact Flash memory storage	Yes	Yes	Yes	No
Power source	4AA batteries, single-use alkaline or NiMH rechargeable	The same as FlexComp Infiniti	The same as FlexComp Infiniti	1 AA battery, single-use alkaline or NiMH rechargeable
Battery life (alkaline)	30 hours typical, 20 hours minimum	The same as FlexComp Infiniti	The same as FlexComp Infiniti	10 hours typical