

Dive into Aquatic Biofeedback™



AquaSense™ comes in a range of sizes for arms and legs.

AquaSense™

Keeping your sensors dry

AQUATIC REHABILITATION USING SEMG™

AquaSense - The first commercially available aquatic biofeedback rehabilitation solution.

Thought Technology's AquaSense is a watertight enclosure designed specifically for use with our SEMG sensors. With AquaSense you can now combine the benefits of aquatic therapy with muscle re-education. AquaSense helps you improve the quality of client treatment, opening new therapeutic opportunities.

Call us regarding upcoming Aquatic Biofeedback workshops in your area!

Product Features:

- Slips easily and comfortably over the limb.
- Made of rugged surgical latex.
- Nonskid grips help prevent slipping.
- Patented integrated pump provides a watertight seal.
- Comes in a variety of sizes for arms or legs.
- Works with all Thought Technology encoders.



Thought Technology Ltd.

2180 Belgrave Avenue Montreal, QC, H4A 2L8, Canada
Tel: 1-800-361-3651 • 514-489-8251 Fax: 514-489-8255
e-mail: mail@thoughttechnology.com
<http://www.thoughttechnology.com>