

BioGraph Infiniti Computer Requirements

Recommended	Minimum
<p>IBM PC compatible (Intel/Pentium/Celeron family or AMD K6/Athlon/Duron family), laptop, or desktop Windows 2000/XP Professional or Home editions. Processor, P4 or higher, 2.4 GHz or higher 40 gigabytes of available hard disk space (For use with 1 hour video recording and processing) Memory, 512 MB of RAM or more CD ROM or DVD drive SVGA graphic card (1024 x 768) or higher resolution adapter & monitor Sound Blaster AWE 32 (or Sound Blaster Live) sound card & Speakers 1 to 4 USB ports, depending on the desired number of Pro/Flex Infiniti encoders 1 to 4 serial com ports, depending on the desired number of ProComp+ encoders. Mouse Desk or Notebook with two monitor capability. CF Reader (For use with compact flash card only) Webcam 30 frames per second (for video purposes only) MS Word 97 or higher (for printing purposes)</p>	<p>IBM PC compatible (Intel/Pentium/Celeron family or AMD K6/Athlon/Duron family), laptop, or desktop Windows 2000/XP Professional or Home editions P4 600MHz Memory, 256MB of RAM or more 1.5 gigabytes of available hard disk space CD ROM or DVD drive SVGA graphic card (1024 x 768) or higher resolution adapter & monitor 16 Bit Sound Blaster compatible sound card & Speakers Mouse 1 to 4 USB ports, depending on the desired number of Pro/Flex Infiniti encoders 1 to 4 serial com ports, depending on the desired number of ProComp+ encoders. Word 97 or higher (for printing purposes)</p>

BioGraph 2.1 Computer Requirements

Recommended	Minimum
<p>IBM PC compatible, laptop, or desktop Windows'98 2000 ***** Processor, P3 1GHz or higher Memory, 256 MB of RAM or more Hard Drive, 30 GB or more Floppy Drive CD ROM SVGA graphic card (16MB) & monitor Sound Blaster AWE 32 (or Sound Blaster Live) sound card & Speakers 1 free serial port, COM1 or COM2 Mouse ***** for XP, please contact our technical support people prior to using XP: 1-514-489-8251 ext 125 or techsupport@thoughttechnology.com *** For split screens: 2 graphic cards and 2 monitors on desktop *** For ProSB Dual: 2 com ports or 2 computers</p>	<p>IBM PC compatible, laptop, or desktop Windows '98,2000 P2 350MHz Memory, 128 MB of RAM Hard Drive, 10 GB Floppy Drive CD ROM SVGA graphic card (4MB) & monitor Sound Blaster compatible sound card & Speakers 1 free serial port, COM1 or COM2 Mouse</p>

CardioPro 2.0 Computer Requirements

<i>Recommended</i>	<i>Minimum</i>
Computer laptop or desktop Windows '98, 2000, XP Processor, Pentium 4 2.4 Ghz Memory, 512 MB of RAM Hard Drive, 40 GB CD ROM 15" SVGA monitor, Video card with 32 or 64 MB Sound Blaster AWE 32 sound card & Speakers 1 free serial port, COM1 or COM2 Mouse	Computer laptop or desktop Windows '98, 2000, XP Processor, Pentium II 350MHz Memory, 128 MB of RAM Hard drive, 5 GB CD ROM 15" SVGA monitor, video card with 4 MB Sound Blaster compatible sound card & speakers 1 free serial port, COM1 or COM2 Mouse

Continence Software Computer Requirements (for MyoTrac3 and MyoTrac3G devices)

<i>Recommended</i>	<i>Minimum</i>
Computer Desktop Windows'98,2000 Processor, Pentium 4 2.4 Ghz Memory, 512 MB of RAM Hard Drive, 40GB available disk space Floppy Drive 17" SVGA monitor, Video card with 32 or 64 MB Sound Blaster AWE 16bit sound card & speakers 1 free serial port, COM1 or COM2 Mouse	Computer laptop or desktop Windows 95, 98,2000 Processor, P2 350 Memory, 128 MB of RAM Hard drive, 5MB available disk space Floppy Drive 15" SVGA monitor, video card with 1 MB Sound Blaster compatible sound card & speakers 1 free serial port, COM1 or COM2 Mouse

EZ Connect Software Computer Requirements (For EEGZ sensor)

<i>Recommended</i>	<i>Minimum</i>
Computer Desktop or Laptop Windows'98, 2000, XP Processor, Pentium 4, 2.4 Ghz Memory, 512MB of RAM Hard Drive, 40 GB disk space available CD ROM 15" SVGA monitor, Video card with 32 or 64 MB Sound Blaster AWE 32 sound card & speakers 1 extra free serial port, COM1 or COM2 Mouse Screen settings at 800x600 pixels	Computer Desktop or Laptop Windows'98, 2000, XP Processor, Pentium 2, 350 MHz Memory, 128 MB of RAM Hard Drive, 10 MB disk space available CD ROM 15" SVGA monitor, Video card with 4 MB Sound Blaster AWE 32 sound card & speakers 1 extra free serial port, COM1 or COM2 Mouse Screen settings at 800x600 pixels