

ProComp 2™

Thought
Technology
Montreal QC Canada

A: EEG ELECTRODE LEAD ONLY
B: EEG, EKG, ACTIVE SENSOR
C&D: EMG, BVP, TMP, SC, RESP, V



EXPAND your
clinical practice!

The Software You Know The Hardware You Trust

The ProComp 2™ is a compact, yet powerful 2 channel version of ProComp+™. ProComp 2™ contains a built-in EEG sensor – requiring only an extender cable for EEG monitoring and biofeedback - and it can use all available ProComp+™ sensors.

ProComp 2™ is available exclusively to owners of ProComp+™ systems and works with your existing BioGraph-registered software.

** Limited Time Offer*



Actual Size

Features

- Monitors 2 channels simultaneously:
 - 2 EEGs
 - 1 EEG and 1 EKG
 - 1 EEG and 1 peripheral
 - 2 peripherals
- Full 13 bit resolution with fiber-optic cable for low noise and flexibility
- Easily worn on headband or shirt collar

Benefits

- Works with your full functioning BioGraph® software
- Uses most of your existing custom settings, screens and protocols
- Quickly set up additional training rooms
- Quality and support from Thought Technology Ltd.

Specifications

- Size: 2.2 x 2.8 x .75" (75 x 55 x 19 mm)
- Battery: 1.5V AA (10+ hours)
- Sensors: EEG, EMG, EKG, SC, TEMP, BVP, Respiration and Voltage.

*If you don't own BioGraph®
Contact us at 1-800-361-3651*



Thought Technology Ltd.

2180 Belgrave Avenue Montreal, QC, H4A 2L8, Canada
Tel: 1-800-361-3651 • 514-489-8251 Fax: 514-489-8255

e-mail: mail@thoughttechnology.com • <http://www.thoughttechnology.com>

MAR 598-01